



creos™



nobelbiocare.com

YOU ARE WHAT
WE STAND FOR.

creos

advanced

Upwards and on the rise

**Innovation is
our passion
but it's your reality
THAT DRIVES US.**

Enhancing the experience

easy-to-use

creos™ is the Nobel Biocare regenerative product portfolio, built to meet your everyday needs. Join us and experience the carefully designed, easy-to-use, effective solutions.

We aim to be the trusted regenerative partner for you and your patients, because you are what we stand for.

Sideways to your side

effective

creos™ xenoprotect

A membrane with outstanding handling that facilitates bone gain



OUTSTANDING HANDLING^{1,2}

- Does not stick to instruments.
- Repositioning in-situ possible.
- Low surface expansion when hydrated.
- Both sides can face the defect.

HIGH MECHANICAL STRENGTH^{2,3,4}


- Excellent suture retention.
- Highly tear-resistant.

ENDURING RESISTANCE TO DEGRADATION IN VIVO^{*3}

- Manufacturing process intended to preserve the natural structure of the collagen fiber network, to confer a high resistance to degradation.

FACILITATES BONE GAIN^{2,3,5,6,7,8}

- Excellent tissue compatibility.³
- Good clinical results.⁵



"What I like is that the handling is very easy. The mechanical stability is very high and when it is rehydrated it adapts very well to the underlying bone"

Dr. Bastian Wessing, Germany

*As shown in an animal model (rat, subcutaneous)




High mechanical strength

In an in vitro study aiming to compare the mechanical strength of commonly used native non-chemically cross-linked and chemically cross-linked collagen membranes⁴:

- creos™ xenoprotect demonstrated the highest force at break, wet (21.2 N).
- creos™ xenoprotect had the highest suture retention, wet (6.1 N).

Comparison of commercial membrane in a hydrated state




Non cross-linked collagen membranes (NXL) – CX: creos™ xenoprotect [Nobel Biocare]; CO: Copios [Zimmer]; JS: Jason [botiss]; OF: Osseoguard Flex [3i]; BG: Bio-Gide [Geistlich]
 Cross-linked collagen membranes (XL) – BE: BioMend Extend [Zimmer]; ML: Mem-Lok [BioHorizons]; OP: OssixPlus [Datum Dental]; BM: BioMend [Zimmer]; CY: Cytoplast RTM [Osteogenics]
 *Statistically significant

Enduring resistance to degradation in vivo without chemical cross linking³


Manufacturing process intended to preserve the natural structure of the collagen fiber network, to confer a high resistance to degradation.³

In an animal model, after 20 weeks, the thickness of xenoprotect decreased only slightly, whereas the reference membrane showed a thickness loss of around 50%, confirming the higher stability of xenoprotect against biodegradation in vivo.³



Representative histological images at 20 weeks implantation in a rat model.

Facilitates new bone formation^{2,3,5,6,7,8}




*Statistically significant

In a comparative in vivo study, creos™ xenoprotect demonstrated significantly higher new bone formation in the central portion of the defect.

This increase in bone formation was associated with significantly increased expression of the growth factor *Bmp2*, which has a strong role in osteogenesis.⁷

In a randomized controlled clinical trial, 24 patients were treated with creos™ xenoprotect and 25 with a reference membrane. In the creos™ xenoprotect group, the defect height reduced at 6-month re-entry by 81%. In the reference membrane group, the defect height reduced at 6-month re-entry by 62%.⁵

Schematic showing the defect height prior to treatment and 6 months after GBR



Scan the code for clinical resources.

creos™ xenogain


3 methods of application to meet all your bone grafting needs



Regenerating bone for 15 years



THREE DIFFERENT METHODS OF APPLICATION:



vial

bowl

syringe

SIMILAR TO HUMAN BONE

- Chemical composition: Ca/P ratio.
- Interconnected macropores.^{9, 10}

EASY HANDLING

- Homogenous particle size.⁹
- Hydrophilic for fast rehydration.^{11, 12}

SOLID FOUNDATION FOR DENTAL IMPLANT TREATMENT

- Osteoconductive properties.¹⁰
- Long-term volume stability.¹³
- Uneventful healing.^{8, 11, 12, 13, 14}

"I appreciated its handling properties and I see its high hydrophilicity as a biological advantage in sinus grafting and peri-implant defect regeneration"

Dr. Werner Zechner, Austria



Bovine

creos™ xenogain collagen



block

syringe

Purified cancellous bovine bone mineral granules and 10% porcine collagen in block form and syringe.

The collagen helps to hold creos™ xenogain collagen in the desired place.

Especially recommended for extraction socket management.



Bovine and Porcine



Scaffold for successful regeneration

Preserved natural features of bone through optimized manufacturing process.¹⁰

CHEMICAL COMPOSITION

With a calcium phosphate ratio that reflects the composition in human bone and a structure with low crystallinity. The body accepts creos™ xenogain as a suitable framework for bone formation.⁹

PARTICLE SIZE



Homogenous particle size.⁹
Maintains space for bone regeneration.¹²

PRESERVED NANOSTRUCTURE

Nanostructure preserved thanks to treatment at comparatively low temperature (600°C) and no sintering.¹⁰

MACRO AND MICRO-STRUCTURE

Interconnected macropores allow cells to invade bone grafts and micropores contribute to capillary liquid uptake (hydrophilicity).^{15, 16}




Photographic micrograph of creos™ xenogain and reference product showing the particle size distribution (magnification 20x)

Solid foundation for implant placement

The graft integrates with the newly formed bone, building a basis for successful implant placement.¹²

Initial situation before GBR 8 months post-surgery




Schematic showing the defect and bone size prior to and after GBR

In a multicenter clinical study involving 46 patients, bone increase after 8 months was 4.0 mm (+56.9 % gain) and 4.7 mm (51.0 % gain) at 1 and 3 mm from the top of the crest, respectively.⁸

GBR led to robust bone regeneration during the 8 months of healing, enabling successful placement of 91 implants in 43 patients, with an average insertion torque of 37.8 ± 5.1 Ncm.⁸

Histological assessment of the trephine cores showed 37.3 % of new bone, 39.1 % of graft material, and 23.6 % of soft tissue (n = 6 cores, 3 patients).⁸




Histological cross section of the cellular components; NB – new bone, XG – graft, scale bar shown in the bottom right corner, red dashed line: bone to graft particle contact.



Scan the code for clinical resources.

creos™ mucogain

Unique oriented porous structure designed for guided soft tissue regeneration



SUBSTITUTES THE NEED FOR A SECOND SURGICAL SITE ^{17, 18, 19}

PATENTED MANUFACTURING METHOD

- Open interconnecting porous structure.
- Designed to promote soft tissue regeneration through the migration of cells and blood vessels into the matrix.^{20, 21, 22}

VARIETY OF CHOICES

- A choice of different sizes and thicknesses.

EXCELLENT HANDLING

- Easy to use.²³
- High suture retention and stress resistance.²³
- Memory effect after hydration and cycling loading in vitro.²⁰
- Trim to precisely fit surgical site.²³

CLINICALLY EFFECTIVE

- Shown to promote soft tissue health and maintain adequate soft tissue thickness in a clinical study.^{23, 24, 25, 26}

"It felt like an autogenous tissue graft and the mechanical stability is amazing"

Dr. Miguel González Menéndez, Spain




Porcine

Use straight out of the box

creos™ mucogain is intended to be used for soft tissue augmentation indications in the oral cavity around teeth or implants:

- Guided tissue regeneration (GTR) procedures in recession defects for root coverage.
- Localized gingival augmentation to increase keratinized tissue around teeth and implants.




Unique oriented porous structure




1. Matrix structure

Interconnecting porous structure, produced by a patented process.^{20, 21, 22}



2. Mechanical properties

After hydration and compression in 49 cycles in vitro, the graft regains its initial volume.²⁰



3. Biological outcome

Designed to promote soft tissue regeneration through the migration of cells and blood vessels into the matrix.^{20, 22}

Clinically effective ^{23, 24, 25, 26}

Clinically effective for soft tissue regeneration in combination with immediate implant placement and bone grafting procedure.^{23, 24}

A retrospective analysis including 45 patients with a follow-up of up to 4.5 years (mean of 1.8 ± 1.3 years) demonstrated that creos™ mucogain promotes soft tissue health and maintains adequate soft tissue thickness when used simultaneously with implant placement.²⁵

Clinical case

Cirillo F. (March 2020).

Periodontal plastic surgery: gingival recession coverage with a xenogenic collagen matrix. The Foundation for Oral Rehabilitation (FOR.org): <https://bit.ly/2TkLsgu>
Images reprinted with permission of the author and FOR.



Buccal view prior surgery (left) and 8 months after surgery (right) on #22, #24, #25, #26 after treatment with creos™ mucogain.



Scan the code for clinical resources.

References

1. Arrighi I, Wessing B, Rieben A, et al. Resorbable collagen membranes expansion in vitro. J Dent Res 2014;93 (Spec Iss B):631
2. Wessing B, Emmerich M, Bozkurt A. Horizontal ridge augmentation with a novel resorbable collagen membrane: a retrospective analysis of 36 consecutive patients. Int J Periodontics Restorative Dent 2016;36(2):179–187.
3. Bozkurt A, Apel C, Sellhaus B, et al. Differences in degradation behavior of two non-cross-linked collagen barrier membranes: an in vitro and in vivo study. Clin Oral Implants Res 2014; 25(12):1403-1411
4. Gasser A, Wessing B, Eummelen L, et al. Mechanical stability of collagen membranes: an in vitro study. J Dent Res 2016;95(Spec Iss A): 1683
5. Wessing B, Urban I, Montero E, et al. A multicenter randomized controlled clinical trial using a new resorbable non-cross- linked collagen membrane for guided bone regeneration at dehiscd single implant sites: interim results of a bone augmentation procedure. Clin Oral Implants Res 2017;28(11):e218–e226.
6. Sanz-Sanchez I, Wessing B, Polizzi G, et al. Randomized clinical trial comparing two resorbable collagen membranes demonstrates good bone formation and soft tissue healing with GBR at single implant sites with dehiscence defects. J Clin Periodontol 2018;45(S19):19–20 [oral presentation].
7. Omar O, Dahlin A, Gasser A, et al. Tissue dynamics and regenerative outcome in two resorbable non-cross-linked collagen membranes for guided bone regeneration: A preclinical molecular and histological study in vivo. Clin Oral Implants Res 2018;29(1):7–19
8. Aleksic Z, Milikovic I, Laziv Z, et al. A multicenter clinical investigation demonstrates bone regeneration in severe horizontal defects in the posterior mandible using creos™ xenoprotect: Interim results. J Clin Periodontol 2018;45(S19):306
9. Nobel Biocare, data on file.
10. Rhee S-H, Park HN, Seol Y-J et al. Effect of heat-treatment temperature on the osteoconductivity of the apatite derived from bovine bone. 2006 Key Engineering Materials 309-311:41-44
11. Park HN, Han, SH, Kim KW, et al. A study on the safety and efficacy of bovine bone-derived bone graft material (OCS-B). J Korean Acad Periodontol. 2005 Jun;35(2):335 – 43
12. Park JB, et al. Maxillary sinus floor augmentation using deproteinized bovine bone-derived bone graft material (OCS-B). Clinical and histologic findings in humans. The Journal of the Korean Dental Association. 2007;45(8):491–99
13. Shin S-Y, et al. Long-term results of new deproteinized bovine bone material in a maxillary sinus graft procedure. J Periodontal Implant Sci. 2014;44:259-64.

14. Kim Y-T, et.al. Periodontal Repair on Intrabony Defects treated with Anorganic Bovine-derived Xenograft.J Korean Acad Periodontol. 2007;37(3):489 – 96
15. Fernandez de Grado G, Keller L, Idoux-Gillet Y et al. Bone substitutes: a review of their characteristics, clinical use, and perspectives for large bone defects management. Journal of Tissue Engineering Volume 9: 1–18, 2018

16. Zhang K , Fan Y, Dunne N et al. Effect of microporosity on scaffolds for bone tissue engineering. Regenerative Biomaterials, 2018, 115–124

17. Aguirre-Zorzano LA, García-De La Fuente AM, Estefanía-Fresco R, et al. Complications of harvesting a connective tissue graft from the palate. A retrospective study and description of a new technique. J Clin Exp Dent 2017;9(12):e1439-e1445.

18. Griffin TJ, Cheung WS, Zavras AI, et al. Postoperative complications following gingival augmentation procedures. J Periodontol 2006;77(12):2070-2079.

19. Harris RJ, Miller R, Miller LH, et al. Complications with surgical procedures utilizing connective tissue grafts: a follow-up of 500 consecutively treated cases. Int J Periodontics Restorative Dent 2005;25(5):449-459.C513

20. Damink L.O., Heschel I, Leemhuis H. et al Gasser A, Wessing B, Eummelen L, et al. Soft tissue volume augmentation in the oral cavity with a collagen-based 3D matrix with orientated open pore structure. Current Directions in Biomedical Engineering 2018; 4(1): 237 – 241

21. Heschel I, et al. 2002. Method for producing porous structures. US patent 6,447,701 B1.

22. Boekema B., Vlig F, Damink O.L. et al. Effect of pore size and cross-linking of a novel collagen-elastin dermal substitute on wound healing. J Mater Sci Mater Med 2014;25(2):423-433

23. Wessing B, Vasilic N. Soft tissue augmentation with a new regenerative collagen 3-D matrix with oriented open pores as a potential alternative to autologous connective tissue grafts. Clin Oral Implants Res. 2014 Sep;25(s10)

24. Sanz-Martin I, Encalada C, Sam-Sanchez I, et al. Soft tissue augmentation at immediate implants using a novel xenogeneic collagen matrix in conjunction with immediate provisional restorations: A prospective case series. Clin Implant Dent Relat Res 2019;21(1):145-153.

25. Wessing B, Boekema B and Gonzalez-Martin. Soft tissue augmentation with a collagen-based 3D matrix with directed pore channels. Clinical Oral Implant Research 2019:Volume30, IssueS19.

26. Cirillo I and Encalada C. Periodontal plastic surgery: treatment of multiple gingival recessions. FOR.org

Products

creos™ xenoprotect

porcine collagen membrane

Size	Article No.
15x20 mm	N1520
25x30 mm	N2530
30x40mm	N3040

creos™ mucogain

porcine collagen matrix

Size	Thickness	Article No.
15x20 mm	3 mm	MU15203
25x30 mm	3 mm	MU25303
15x20 mm	5 mm	MU15205
25x30 mm	5 mm	MU25305

creos™ xenogain

deproteinized bovine bone matrix

Weight	Granule size	Volume	Vial	Bowl	Syringe
0.25 g	Small (0.2-1.0 mm)	0.36 cc	N1110	N1110-B	N1210
	Large (1.0- 2.0 mm)	0.54 cc	N1111	N1111-B	N1211
0.50 g	Small (0.2-1.0 mm)	0.82 cc	N1120	N1120-B	N1220
	Large (1.0- 2.0 mm)	1.27 cc	N1121	N1121-B	N1221
1.00 g	Small (0.2-1.0 mm)	1.71 cc	N1130	N1130-B	
	Large (1.0- 2.0 mm)	2.69 cc	N1131	N1131-B	
2.00 g	Small (0.2-1.0 mm)	3.64 cc	N1140	N1140-B	
	Large (1.0- 2.0 mm)	5.74 cc	N1141	N1141-B	

creos™ xenogain collagen

creos™ xenogain + 10% porcine collagen type I

Size	Block size	Article No.
100 mg	6 x 6 x 6 mm	N1320
250 mg	7 x 8 x 9 mm	N1330
500 mg	9 x 10 x 11 mm	N1340

Size	Syringe size	Article No.
250 mg	4.6 x 40 mm	N1410
500 mg	5.6 x 45 mm	N1420

YOU ARE WHAT
WE STAND FOR.

ORDER ONLINE

Order the complete creos™ range
24 hours a day through the
Nobel Biocare online store.

nobelbiocare.com/store

ORDER BY PHONE

Call our customer service team, or
contact your sales representative.



nobelbiocare.com



87012 GB 2012 Printed in the EU © Nobel Biocare Services AG, 2020. All rights reserved. Distributed by: Nobel Biocare. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information.

Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Please contact the local Nobel Biocare sales office for current product assortment and availability.

Disclaimer: For prescription use only. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.

Legal manufacturer for creos™ xenoprotect and creos™ mucogain is Matricel GmbH, Kaiserstrasse 100, D-52134 Herzogenrath, Germany. Legal manufacturer for creos™ xenogain and creos™ xenogain collagen is NIBEC Co., Ltd., Iwol electricity-electronic Agro-Industrial Complex, 116, Bamdi-gil, Iwol-myeon, Jincheon-gun, Chungcheongbuk-do, 27816, Korea.